


PETTINARI



FLUOROSCOPIC GENERAL SURGERY TABLE MILLENIUM®
APT FOR "C ARC" IN ALL THE ESPECIALITIES,
WITH MIXED SYSTEM WHICH LEAVES THE NON TRANSPARENT
SECTION OUT OF WORK.
CHANGEABLE HEIGHT BY PEDAL HIDRAULIC SYSTEM
WITH STAINLESS STEEL TELESCOPIC COVER
TRENDELENBURG AND REVERSE + 25° - 15°
LEFT AND RIGHT LATERAL TILT + 15° - 15°
BOTH CONTROLLED BY BILATERAL MECHANICAL SYSTEM



FIVE SECTIONS RADIOLUCET TABLE TOP
BACK / ABDOMINAL AND PODALIC SECTION
ARTICULATED BY PNEUMATIC SYSTEM
FULL STAINLESS STEEL ESTRUCTURE
FULL STAINLESS STEEL BASE
WITH FIXING SYSTEM FOR RETRACTILE
WHEELS BY PEDAL
ACCESORIES INCLUDED:
(1) ARTICULATED NARCOSIS FRAME
(1) PAIR OF REGULABLE SHOULDER HOLDERS
(1) PAIR OF REGULABLE SIDE SUPPORTS
(1) SET OF FOAM MATTRESS WITH PLASTIC SHEET
MOD.: PM 220 HNA1

