

PETTINARI



INTERNMENT BED PM NOVA®

STEEL SHEET LITTER SAE 1010
ANTIRUST TREATMENT
FIXED HEIGHT

24 mm. (1") TWO SIDED "MASISA" SHEET HEAD AND
FOOT - BOARD, ABS EDGES

(2) HARD RUBBER BUMPERS

NEW POSITIONS - VASCULAR AND FOWLER POSITIONS
MECH - MANUALLY OPERATED

(4) ADJUSTABLE STAINLESS STOPPERS

MEASUREMENTS: 1,95 X 0,83 X 0,62 m (6,5 X 2,8 X 2 ft)

MOD.: PM 100 HNA1 - ADULTS -

MOD.: PM 120 HNA1 - PEDIATRICS -

INTERNMENT BED PM CENTURY®

STEEL SHEET LITTER SAE 1010
ANTIRUST TREATMENT
FIXED HEIGHT

24 mm. (1") TWO SIDED "MASISA" SHEET HEAD AND
FOOT - BOARD, ABS EDGES

(2) HARD RUBBER BUMPERS

NEW POSITIONS - VASCULAR AND FOWLER POSITIONS
MECH - MANUALLY OPERATED

(4) ADJUSTABLE STAINLESS STOPPERS

MEASUREMENTS: 1,95 X 0,83 X 0,62 m (6,5 X 2,8 X 2 ft)

MOD.: PM 104 HNA1 - ADULTS -

MOD.: PM 124 HNA1 - PEDIATRICS -



INTERNMENT BED PM MILENIUM®

STRATIFIED LITTER ANTIDECAJ SINTRA SHEET

PROFILE FRAME MILENIUM DGA

NEW POSITIONS - VASCULAR AND FOWLER POSITIONS

45°TOTAL TRENDELENBURG AND INVERSE

ELECTRONICALLY OPERATED

REMOTE CONTROL

220 VCA - 50 Hz IMPUT

HEAD AND FOOT - BOARD MADE IN "MASISA" SHEET ON
BOTH SIDES AND FINISHED IN ABS EDGES, TUBULAR

PERIMETRAL FRAME SAE 1010 STEEL, REINFORCED WITH

(4) CIRCULAR RUBBER BUMPERS

(2) MANUALLY FIXED FOLDING LONG LATERAL BED - RAILS
AUTOMATIC HOOK UP

(1) STAND WITH (3) MULTIPURPOSE ADJUSTABLE HOOKS
TO HANG DRIPS, BLOODBAGS, ETC.

(4) 100 mm. (4") DIAMETER DIRECTABLE WHEELS TOTAL BLOCK

MEASUREMENTS: 1,95 X 0,83 m (6,5 X 2,8 ft)

0,50 m (1,7 ft) MIN. HEIGHT

0,75 m (2,5 ft) MAX. HEIGHT

MOD.: PM 126 HNA1



ROLLING UNIT FOR KIDS PM NOVA®

RIGID METALIC WIRENING WITH FOWLER ARTICULATION

MECH - MANUALLY BACKREST

(2) RUBBER BUMPERS

(1) PANORAMIC FOOT - BOARD

(4) SLIDING LATERAL PANORAMIC BED - RAILS

(4) U - SHAPED TUBULAR LEGS

(4) 100 mm (4") DIAMETER DIRECTABLE WHEELS WITH TOTAL
BLOCK

MEASUREMENTS: 1,30 X 0,60 X 0,70 m (4,3 X 2 X 2,3 ft)